


THREE NECK BOTTOM OUTLET SPHERICAL VESSELS



Cat. Ref.	Nominal Capacity	L	DN	DN1	DN2
VSM5	5 L	425	50	25	25
VSM10	10 L	575	80	25	25
VSM20*	20 L	650	100	25	25
VSM50*	50 L	840	150	40	40
VSM100	100 L	950	225	40	40
VSM200	200 L	1200	300	40	40
VSM300	300 L	1310	400	50	50
VSM500	500 L	1450	450	50	50