


SPHERICAL VESSELS WITH WIDE BOTTOM OUTLET



These vessels are generally used to fit immersion exchangers in the bottom. Special heating mantle or bath should be used if used with.

Cat. Ref.	Nominal Capacity	L	DN	DN1	DN2	DN3
VSR50	50 L	790	150	150	40	100
VSR100	100 L	900	225	150	40	100
VSR200	200 L	1150	300	150	40	100
VSE50	50 L	840	150	225	40	100
VSE100	100 L	950	225	225	40	100
VSE200	200 L	1200	300	225	40	100

* marked items are available fast.